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INSTRUCTIONS

Case content:

- CANAL SEALER SYSTEM POWDER (8g): Zinc Oxide, Titanium Oxide, Zinc Stearate, Calcium Wolframate, Zinc Acetate, Preservative Complex “Obcol 5”.
- CANAL SEALER SYSTEM LIQUID “PLASTIC” (Eugenol-free) (5ml): Rosin, EBA, Biosol.
- CANAL SEALER SYSTEM CREAM (1ml): Surfactant, Preservatives, PEG, Distilled Water.
- CANAL SEALER SYSTEM WASH (100ml): Detergent hydro-alcoholic solution. Preservatives.
- CANAL SEALER SYSTEM TCF (0.5ml): Preservative Complex in not water-soluble cream (Vaseline – lanoline).

CANAL SEALER SYSTEM (C.S.S.)

Definition:

C.S.S. CEMENT is a permanent endodontic, root canal cement for use after vital pulpectomy or after necrotic, gangrenous canal treatment. *Only for dental professional use.*

Advantages:

This filling creates the best conditions for subsequent periapical osteo-cementum formation (so-called biological “retrograde” seal).

C.S.S. Powder mixed with C.S.S. Liquid “Plastic” assumes, after setting, a gutta-percha similar consistency, substituting in fact all its physical properties. In addition it has the advantage of a preservative effect, greater adhesion, better saturation and no expansion.

Indications:

- The best results can be obtained with the sequential application of the methodology: C.S.S. Cream - C.S.S. Wash – C.S.S. TCF Temporary (if necessary), followed by the final, permanent filling with C.S.S. Cement. Nevertheless other traditional sealing methods can also be used before the final filling with C.S.S. CEMENT.
- The cement is radio-opaque.

Contraindications:

- Avoid any prolonged contact with the mucous membrane by protecting the site appropriately (rubber-dam, cotton rolls and, if necessary, rinse with water).
- In case of eugenol allergy use definitely the “Eugenol-free” C.S.S. Liquid Plastic.
- In case of composite or Glass-ionomer crown filling use certainly “Eugenol-free” C.S.S. Liquid Plastic.

Cautions:

- The final permanent canal filling is correct when it reaches the apical constriction (Foramen physiologicum).
- In spite of the high bio-compatibility, it’s appropriate to avoid “bringing” the C.S.S. products beyond the anatomic foramen (apex).
- During the canal obturation insulate the “site” with cotton rolls or rubber-dam.
- If permanent filling has to be delayed, the canal should be protected with an intermediary preservative cement (See: TCF- Temporary Canal Filling).
- In open apices it is suggested to delay permanent seal and place a temporary prolonged protective filling (awaiting the conclusion of natural apexification).
- A temporary protective filling is also indicated in the presence of risky zones (mandibular canal, maxillary sinus) where inadvertent over fillings must be avoided. The final permanent filling follows only after biological apical seal, performed by osteo-cementum.

Mixing Instructions:

- Put **one** (1) drop of C.S.S. Liquid Plastic on a sterile glass slab (or mixing pad) and add C.S.S. Powder progressively using a flexible spatula. Spatulate to obtain a creamy mix.
- Working time is about 5 to 6 minutes.

Notice:

- Dispose of any waste as special, not dangerous, garbage.
- Keep away from children.
- Use sterile tools.
- Special precautions: None.
- Undesirable secondary effects: None at the current state of knowledge.

- Cloth cleansing: remove mixed cement only after complete setting.

Storage: Store at room temperature. No need for refrigeration.

OPERATIVE-AID PRODUCTS IN SEQUENTIAL PROGRAM

C.S.S. Cream: Lubricant, Surfactant, Preservative Effect, Water-soluble cream.

Indications: Root canal preparation with files and reamers (Lubricating effect in normal and narrow canals).

Debris and Smear-layer removal.

Advantages: Ready for use with preservative effect for short-term application. The lubricating effect facilitates root canal instrumentation. The cream, in contact with dentinal tubulis and secondary canals, functions as protecting preservative medium. It englobes debris and smear-layer, inhibiting bacterial presence and facilitating their removal from canal.

Directions: Take a minimal dose from dispensing syringe and put it in a sterile Dappen-Dish, (glass slab or mixing pad). Coat all canal instruments before use. At the end of instrumentation remove the used cream with final uncoated instrument and rinse with C.S.S. Wash or sterile water.

Caution: The C.S.S. Cream is hydro-soluble and can be easily removed with C.S.S. Wash.

Contraindications: not ascertained.

Storage: Store at room temperature below 30° C. If product is too fluid refrigerate temporarily to get back its original thickness.

C.S.S. Wash: Irrigating detergent Hydro-Alcoholic solution.

Indications: Root canal irrigation ready for use. C.S.S. Cream solvent.

Advantages: C.S.S. Wash has a detergent action. It helps eliminate the cream (containing englobed debris). The presence of alcohol helps to obtain a clean dry field, ready to receive the final filling with C.S.S. Cement.

Directions: Fix disposable canalar tip to pipette. Insert the tip into dispensing bottle and fill pipette. Insert tip into canal orifice and irrigate without pressure.

Caution: C.S.S. Wash contains Alcohol, then it is flammable. Avoid direct exposure to flames.

Contraindications: None at the current state of knowledge.

Storage: Store at room temperature. To avoid evaporation keep bottle well closed.

C.S.S. TCF - Temporary Canal Filling (Prolonged): Not water-soluble cream with preservative effect for prolonged-term application (4 weeks). Re-applyable. For a temporary use.

Indications: - Temporary canal filling in all those cases where permanent obturation has to be delayed.

- Apexification (presence of an incomplete apical closure, open apices).
- Presence of risky zones (proximity of mandibular canal, maxillary sinus) where inadvertent overfilling of endodontic cement must be avoided.
- The final permanent filling follows biological osteo-cementum apical seal.

Advantages:

- Non hardening vehicle. The cream is ready for use. It is stable being not water-soluble.
- Local protective preservative effect. Reapplicable.
- Not radio-opaque: permits periapical healing procedures control (evolution of apexification).
- No odour. Does not colour the tooth.

Directions: After root canal preparation deposit a small amount of cream, inserting the tip into canal. The canal should be completely filled by pushing the syringe plunger while slowly withdrawing the tip. Cream can be pushed into canal by a paper cone as well.

Caution: Do not overfill in periodontium. Do not inject!

Contraindications: None at the current state of knowledge.

Notice: Removable with Ether or Chloroform.

Storage: Store at room temperature. No need for refrigeration.

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